








Tav. N°
C

SEZ.A:A



PART. X
RETE IN PA 6.6
SCALA 10:1



ORIGINAL IN RED

NOTA: SUPERFICIE FILTRANTE 5.52 cmq

DIMENSIONAL TOLERANCES
ACCORDING TO EN 20286
(ISO 286)

SCALE
SCALA

2:1

EN 20286	Js 01	Js 4	Js 9	Js 14
	Js 0	Js 5	Js 10	Js 15
	Js 1	Js 6	Js 11	Js 16
	Js 2	Js 7	Js 12	Js 17
	Js 3	Js 8	Js 13	

DATE
DATA

25/01/2008

CHECKED - VISTO *S.* *A.*
CHANGED FROM - MODIFICHE

	(2) RETE	7410202700	MEDIFAB 263 μm PA-6.6 MONOFILAMENT		
	(1) SUPPORTO	9700650000	TERLUX 2802 HD ABS		
	12408500				
ITEM	ARTICLE CODE	MATERIAL CODE	TITLE	MATERIAL OR TITLE	NOTE
			FILTRO SECURA 263 μm		DRAWING CODE NR. STAMPO 124

THIS MATERIAL IS THE PROPERTY OF LUCOVED S.p.A. THE INFORMATION IS CONFIDENTIAL AND IS TO BE USED ONLY IN CONNECTION WITH MATTERS AUTHORIZED BY LUCOVED S.p.A. AND NO PART OF IT IS TO BE DISCLOSED TO OTHERS WITHOUT PRIOR WRITTEN PERMISSION FROM LUCOVED S.p.A.