




THIS MATERIAL IS THE PROPERTY OF LUCOMED S.p.A. THE INFORMATION IS CONFIDENTIAL AND IS TO BE USED ONLY IN CONNECTION WITH MATTERS AUTHORIZED BY LUCOMED S.p.A. AND NO PART OF IT IS TO BE DISCLOSED TO OTHERS WITHOUT PRIOR WRITTEN PERMISSION FROM LUCOMED S.p.A.

Tav. N°
C

ZIGRINATURA
N°30 DENTI




SEZ. A:A




FILETTATURA $\frac{1}{2}$
PRINCIPI $\frac{1}{2}$ PASSO

$\varnothing 7.00 \pm 0.10$
 $\varnothing 7.90 \pm 0.10$



**ORIGINAL IN
RED**



CONFORME ALLE NORME ISO-594

DIMENSIONAL TOLERANCES
ACCORDING TO EN 20286
(ISO 286)

SCALE
SCALA

4:1

33628503L

A 9700650000
B 9700600300

TERLUX
2802 HD
SINKRAL
L 322

ABS
ABS

+2% Col. BIANCO
Cod. 9780030000
17178 BIANCO

EN 20286	Js 01	Js 4	Js 9	Js 14
	Js 0	Js 5	Js 10	Js 15
	Js 1	Js 6	Js 11	Js 16
	Js 2	Js 7	Js 12	Js 17
	Js 3	Js 8	Js 13	

DATE
DATA
25/01/2008

CHECKED - VISTO
CHANGED FROM - MODIFICHE

ARTICLE DESCRIPTION
DESCRIZIONE
**TAPPO ABS BIANCO
LUER LOCK MASCHIO**

DRAWING CODE
NR. STAMPO

**336-2L
ABS-W**