


THIS MATERIAL IS THE PROPERTY OF LUCOMED S.p.A. THE INFORMATION IS CONFIDENTIAL AND IS TO BE USED ONLY IN CONNECTION WITH MATTERS AUTHORIZED BY LUCOMED S.p.A. AND NO PART OF IT IS TO BE DISCLOSED TO OTHERS WITHOUT PRIOR WRITTEN PERMISSION FROM LUCOMED S.p.A.






Tav. N°
C

ZIGRINATURA
N°30 DENTI



**ORIGINAL IN
RED**



CONFORME ALLA NORMA ISO-594

DIMENSIONAL TOLERANCES
ACCORDING TO EN 20286
(ISO 286)

SCALE
SCALA

4:1

EN 20286	Js 01	Js 4	Js 9	Js 14
	Js 0	Js 5	Js 10	Js 15
	Js 1	Js 6	Js 11	Js 16
	Js 2	Js 7	Js 12	Js 17
	Js 3	Js 8	Js 13	

DATE
DATA
25/01/2008

CHECKED - VISTO *[Signature]*
CHANGED FROM - MODIFICHE *[Signature]*

ARTICLE DESCRIPTION
DESCRIZIONE
**RACCORDO $\varnothing 2.2$
LUER LOCK MASCHIO**

DRAWING CODE
NR. STAMPO
455L

ITEM	ARTICLE CODE	MATERIAL CODE	TITLE	MATERIAL OR TITLE	NOTE
	45508500L	9700650000	TERLUX 2802 HD	ABS	